

TEST FIXTURE TO DETERMINE FORCES TO OPERATE PISTON

of sterile hypodermic syringes for single manual use according to DIN EN ISO7886-1 Fig.E1


- consists of holder and inserts suitable for testing syringes of different diameters
- incl. height adjustable cup table

Item no.:	THS727-80
Max. load:	50 N
Coupling:	Af159 or Af20. Further couplings on request.
Holder and inserts:	Stainless steel
Cup table:	Aluminium, anodized
Temperature range:	0 ... +70°C Further temperature ranges on request
Weight:	2.7 kg per holder incl. insert and cup table
Scope of supply:	1 fixture



Example

THS727-80 holder with inserts for testing syringes of 3 ml, 10 ml and 20 ml:





ISO 7886-1 fig.E.1

ISO 7886-1 fig.E.1

ISO 7886-1 fig.E.1






Example for special version:



THS727-80-H500

Holder height 500 mm

Item no.:	THS727-100
Adjustable opening:	4 – 38 mm
Coupling:	Af159 or Af20. Further couplings on request.
Fixture material:	Stainless steel
Clamps:	Rubber-coated
Cup table:	Aluminium, anodized
Temperature range:	0 ... +70°C Further temperature ranges on request
Weight:	2.4 kg per holder incl. cup table
Scope of supply:	1 fixture. Cup table (item no. <i>THS727-100-Tbl</i>) is to be ordered separately.



See [Animation](#)